


Spring Best in Class	Company		Animas	Roche Accu Chek	Roche Accu Chek2	Medtronic	Medtronic2	Orbit 90	FIFTY50	Cozmo
	Model	Universal™	Inset II™	Ultraflex	FlexLink	Quick-Set®	Mio	Orbit 90	Thinset	Cleo90®
	Connector Type	Luer-Lock	Luer-Lock	Luer-Lock	Luer-Lock	Luer-Lock, Paradigm	Luer-Lock, Paradigm	Luer-Lock	Luer-Lock	Luer-Lock
	Cannula Length	6, 9 mm	6, 9 mm	8, 10 mm	8, 10 mm	6, 9 mm	6, 9 mm	6, 9 mm	6, 9 mm	6, 9 mm
	Gauge	28 gauge intro needle, 22 gauge cannula	27 gauge intro needle, 0.68mm (width) cannula	27 gauge intro needle, 25 gauge cannula	22 gauge cannula	27 gauge intro needle, 22 gauge cannula	27 gauge intro needle, 0.68mm (width) cannula	26 gauge cannula		28 gauge intro needle, 25 gauge cannula
	Tubing Length	24, 31, 43 inches 60, 80, 109 cm	18, 23, 32, 43 inches 46, 60, 80, 109 cm	24, 31, 43 inches 60, 80, 109 cm	12, 24, 31, 43 inches 30, 60, 80, 109 cm	18, 23, 31, 43 inches 46, 60, 80, 109 cm	18, 23, 32 inches 46, 60, 80 cm	18, 24, 30, 42 inches 46, 60, 76, 109 cm	23, 43 inches 60, 109 cm	24, 31, 42 inches 60, 80, 109 cm
	Inserter	all-in-one, retractable needle	all-in-one	inserter needs to be connected to the infusion set	inserter needs to be connected to the infusion set	inserter needs to be connected to the infusion set	all-in-one	none	none	all-in-one
	Disconnect Type	at site, 360° rotation	at site	at site	at site	at site	at site	at site, 360° rotation	at site	at site, 360° rotation
	Spring Universal BEST-in-CLASS Feature Summary	Retractable, hidden needle - safe/user never needs to see or handle; Detach-Disconnect detection; Small all-in-one inserter - one-handed activation with a push of a button; Use with any luer-lock insulin pump reservoir, Portable sterile and compact packaging including both the inserter and the infusion set; Self-adhesive hypoallergenic tape * The information provided here is based on product labeling at the time it was published								

More information: www.springnow.com