

Spring™ Universal Infusion Set

Feature	Specifications	Benefit
Inserter	Insertion velocity 12 meters/sec.	Lightning speed of insertion minimizes sensation
	37mm (H) x 27mm (W) x Ø28mm (D) of base.	Smallest* in class.
	All-in-one audible one-button push insertion.	One-handed activation with the push of one button. Audible confirmation. Stabilizes insertion. No assembly required.
	Protective covering.	Prevents accidental activation.
Intro Needle	Hidden at all times.	Eliminates sharps injuries and biohazard potential.
	28 gauge auto-retractable needle with a spiral trajectory.	Offers one of the least painful insertion experiences.
	Trocar (rather than bevel) tip.	Could reduce potential for tissue irritation and scarring, as well as minimize risk for cannula bending.
	The needle spring is pre-cocked, 2nd spring retracts the needle.	One of the thinnest subcutaneous insertions available
Cannula	90° straight-angle Teflon. 6, 9 mm = (0.23, 0.354inches) length. 22 gauge.	Versatile cannula depth for individual preferences.
	Lubricated with medical silicone fluid.	Reduces friction and thus sensation.
Disconnect Type	At-site disconnection, 360° rotation for insertion options and placement adjustment, locks in place during use.	Best-in-class for easy temporary disconnection and on-site rotation
	Includes protective cap.	User-friendly for insertion site protection after disconnection.
Base	Detach-Detect mechanism.	The pump alarm sounds upon sensing occlusion resulting from a detachment of the infusion set base thus preventing insulin delivery disruption
	Transparent window for site inspection.	Effortless inspection of the insertion site.
	Water-proof, hypoallergenic self-adhesive tape.	User-friendly.
	Base is attached to the body <i>prior</i> to insertion, a “tunnel” guides the needle and cannula.	Ensures accurate 90 degree insertion. Reduces risk of crimping cannula. Minimizes Sensation.
Tubing Length	24, 31, 43 inches = 60, 80, 109 cm.	Convenient for all body builds and user preferences.
	Co-extrusion of two proprietary materials.	The inner lining is made of a material ensuring tubing durability, the outer lining - flexibility.
	High transparency.	Enhances visual air bubble detection.
Connector	Universal connection to any luer-lock insulin pump reservoir.	Enjoy the benefits of Spring Universal Infusion Set with your favorite insulin pump.

* Based on the internal data comparing leading infusion sets on the market.